

ENGINEERING CASE LIBRARY

Development of a Dynamic Seal

at

Beckman Instruments, Inc. - Spinco Division

Mr. Chet Arnessen*, a young graduate of San Jose State College and an employee of Beckman Instruments for the past several years, was asked to redesign part of the "buffer pump" of his company's Model 130 Spectrochrom Analyzer. A few of Beckman's customers had complained that the moving seals in this pump developed leaks after very little use; these complaints were reaching the company at the rate of one or two a month with some repeats. Beckman's customers were scattered all over the world and, since the Spectrochrom Analyzer is not portable, the company had to send its Field Service personnel to them to make whatever repairs were necessary.

Beckman Instruments

Beckman Instruments, Inc., has been manufacturing instruments for research and industrial use since 1935. The Spinco Division, located in Palo Alto, California, is a principal manufacturer and developer of instruments for biochemical and physiological research. Their instruments are used primarily for ultracentrifugation, electrophoresis, ion-exchange chromatography, microcalorimetry, viscometry, and biomedical measuring and monitoring.

All Beckman Spinco Division instruments are covered by a one year warranty which provides for both preventive maintenance and emergency service. The company's Field Service engineers are responsible for these operations. They instruct customers on the use of Beckman instruments, keep records of all repairs made on them, and make recommendations to the home office on improvements that they feel should be made.

* Fictitious Name

(c) 1967 by the Board of Trustees of the Leland Stanford Junior University.

Prepared in the Design Division of the Department of Mechanical Engineering by William J. Clemens and David Horine under the direction of Professor H. O. Fuchs with financial support from the National Science Foundation.

The Spectrochrom Analyzer

The Model 130 Spectrochrom Analyzer, described in Exhibit 1, is primarily a research instrument employed in the analysis of certain biochemical materials. One of the principal components of this instrument is a gradient buffer pump (Exhibit 2a) which mixes two buffer solutions and pumps the resultant mixture to the chromatograph column. The pump itself consists of three cam-driven plungers operating in cylinders where the buffer solutions, usually organic solvents such as acetic acid, ammonium acetate, and pyridine acetate, are mixed and pumped. The proportion of the two solvents in the mixture is regulated by an interchangeable program cam.

Soon after the first models of the Spectrochrom Analyzer were marketed, Beckman Instruments began to receive complaints that the plungers in the buffer pumps were leaking. The first few were repaired routinely in the field. As a pattern of failure became evident, the decision was made to review and redesign the moving seals associated with these plungers. It was at this time that Chet was assigned the problem.

Design Analysis

Chet began his work by systematically reviewing the design history of the buffer pump. What follows is a brief summary of his findings.

The plunger and cylinder used on the first breadboard design had been a glass hypodermic syringe (shown in Exhibit 3, along with the first stainless steel plunger) with an inside diameter of 0.351 inches. As the buffer solution evaporated, crystals deposited on the ground glass cylinder wall and the piston jammed.

Successive models were tested (see Exhibits 4 and 5) using a variety of piston and cylinder materials with o-ring seals. The most effective proved to be a smooth bore (honed to 4-6 min. finish) stainless cylinder with a stainless steel piston equipped with a Rulon¹ tip and a buna "N" o-ring. After evaluation by the Beckman Applications Research Section, this design was released for production.

Redesign

After he had acquainted himself with the design history of the pump, Chet examined the production model and talked with Field Service Department personnel in an attempt to determine the causes of failure. Most of the people he consulted felt that the problem lay in the type of seal chosen for the plunger. They found that the rubber o-ring was badly abraded on all of the pumps they had checked. There were several opinions expressed about the cause of this condition, but Chet thought that two modes of failure seemed most likely: 1) some of the organic solvents might be causing the o-rings to deteriorate; and 2) a thin film of fluid remaining on the cylinder wall might dry and leave an abrasive surface to wear down the o-rings.

¹"Rulon" is a trade name for a filled Teflon material that has higher wear resistance and less creep than pure Teflon.

Chet discovered that there was very little information available on small seals which could be directly applied to his project. Furthermore, he was unable to locate any seals of exactly the size he needed. "When I found there were no standard sized seals I could use," he explained, "I thought of changing the cylinder and plunger to a standard size, but that would have meant changing the stroke of the pumps in production and those already in the field. This would have been expensive, so we decided to retain the existing diameters on our redesign." While he was looking through one of the manufacturers' catalogues, he found that some data on o-ring performance and mounting recommendations were obtained by mathematically scaling the empirically derived data from one size of o-ring to fit a number of smaller ring sizes. From another manufacturer he learned that industrial recommendations for o-rings are usually acceptable for static seals, but they are not always applicable to dynamic sealing problems.

Chet designed a new cylinder and two plunger seals that did not require o-rings. He replaced the stainless steel cylinder with one of glass; unlike the original breadboard component, however, this cylinder had a precision bore (4-6 min. finish) inner surface.

The first of Chet's two plunger designs is shown in Exhibit 6. It consists of a Teflon tip attached to a "Kel-F"² plastic shank. Three "wiper rings" were machined into the Teflon tip to act as seals and guides. By providing an interference fit between the wiper rings and the glass cylinder, he hoped to minimize leakage past the inert Teflon. After three days of continuous testing, using water as the test solution, the wipers began to leak. Chet measured the outside diameters of the wipers and found that they had shrunk. This was apparently due to cold flowing of the Teflon.

Chet's second plunger seal redesign is shown in Exhibit 7. He had a "lip seal" machined into a Teflon tip which had an interference fit with the glass cylinder wall. This seal was intended to expand outward against the cylinder wall when pressure was applied to the chamber fluid. This design, however, began to leak soon after it was installed in the pump. Chet believed that this was due to cold flowing of the Teflon which had caused the lip seal to pull away from the chamber walls at low pressures.

Chet next tried o-rings made of Teflon. Since these rings could not be stretched into mounting grooves as could conventional rubber o-rings, they could not be tested on the existing steel plunger. Chet devised a test plunger (Exhibit 8) which would accommodate several types of Teflon rings. The Kel-F plastic tip of the plunger was threaded so that it could be screwed down against the o-rings to hold them on the plunger. A spiral Teflon wiper ring was mounted around the middle of the plunger to act as an additional guide and seal. Chet chose the spiral Teflon ring because 1) the Teflon would not swell upon contact with a buffer (as does rubber); and 2) a spiral ring could be installed on the shank of the plunger without the necessity of fabricating it to screw apart to accommodate an ordinary o-ring.


²"Kel-F" is a trademark for a relatively inert plastic.

Chet tested three types of rings: a Teflon coated rubber ring, a slotted Teflon ring with a spring insert to provide a sealing force (Exhibits 7b and 7c), and a plain slotted Teflon ring. The Teflon coated rubber o-ring began to swell after two days of continuous operation. After disassembling the plunger, Chet found that the Teflon coating had worn off, thus exposing the base rubber to the buffer solution. The plain slotted o-ring began to cold flow and leaked soon after it was installed. The slotted ring with a spring insert leaked due to insufficient spring strength to affect sealing.

The seal Chet finally decided upon is the Sealpruf³ Kap Seal shown in Exhibit 9. It makes use of a Teflon cap for the sealing surface and an inner synthetic rubber o-ring that forces the Teflon cap against the cylinder wall. Instead of allowing a clearance between the plunger groove and the seal as the manufacturer recommended for normal sealing problems (Exhibit 9), Chet screwed the plunger cap against the ring as shown in Exhibit 10a. In this design, the Teflon is contained by the mounting slot and thereby prevented from cold flowing.

Chet life-tested the Kap Seal for a month before turning it over to production. The final design, shown in Exhibit 10, was subsequently passed by the Applications Research Section and put into final production. After almost one year in service, no serious flaws in the pump have been reported.

³Fictitious name.




Model 130 Spectrochrom® Analyzer

The Spectrochrom Analyzer offers a complete facility for column chromatography in one instrument . . . providing automated analyses of proteins, peptides, nucleic acids and their derivatives, and other biochemical materials. In a single run, with little or no operator supervision, it prints out data that formerly required many time-consuming consecutive runs, and much supervision. Thus the Spectrochrom can simplify chromatographic research programs, and extend and amplify them as well.

With its inter-related combination of components, the Spectrochrom enables the operator to vary selectively the common parameters of a chromatographic analysis until ideal conditions are established. These can then be reproduced faithfully in following runs.

The sample can be studied at wavelengths that best suit its absorption pattern—e.g., nucleic acid derivatives at 250, 260 and 280 $m\mu$. pH and conductivity gradients can be both programmed and monitored, with readings printed out on the recorder chart for correlation with absorption peaks. The recorder chart is keyed by an event marker to the collected fractions, which are held under refrigeration ready for further analysis.



Major Components

Jacketed analytical or preparative columns can be mounted on the Model 130, and cooled or heated from 4° to 75°C. Quick-disconnect fittings permit fast column interchange.

Buffer gradients are supplied by Beckman-designed buffer pump, with gradient programmed by the operator as desired.

Effluent pH and conductivity monitors provide a printed record of these measurements—facilitating the choice of operating conditions.

Dual-beam spectrophotometer automatically scans column effluent at choice of three wavelengths from 220 to 750 $m\mu$; two cuvette path lengths provide excellent sensitivity for a wide range of sample concentrations.

10-Channel recorder provides six color-coded curves of effluent absorbance; two other channels are for pH and conductivity monitoring, leaving two for auxiliary uses.

Refrigerated fraction collector contained in the Spectrochrom Analyzer is the Model 132 Collector described on the following page. It is held in a refrigerated drawer that keeps the collected fractions cold.

Controls permit the programming of a change of column temperature, change of buffer, and instrument shut-down. Panel connections and a valving system allow column effluent to be diverted through an auxiliary instrument and back to the Model 130's analytical system.

318003 Model 130 Spectrochrom Analyzer, 68" x 30" x 76" high, with gradient pump, spectrophotometer, pH and conductivity meters, flowmeter, fraction collector, and controls.

Exhibit 1.




Exhibit 2a. Gradient Buffer Pump

Insert; Rear view of linkage showing push-stop for program cylinder arm.




Exhibit 2b. Buffer Pump Flow Diagram




Exhibit 3. Hypodermic syringe, showing glass cylinder, glass plunger, and first stainless steel plunger.




Exhibit 4. First production model stainless steel cylinder and plunger. Hexagon fitting was added to make cleaning easier. A "Rulon" insert was added as shown in Exhibit 5 to provide a second contact surface for the plunger to ride on.




Exhibit 5. Exploded view drawing of cylinder and plunger shown in Exhibit 4.




Exhibit 6. "Wiper ring" plunger seal design. Rings are machined on the Teflon tip as shown on the sketch below.





Exhibit 7. So-called "lip seal" plunger seal design. A cross-section of the Teflon tip is shown below.






Exhibit 8a. Assembled test plunger. Note the spiral wiper ring installed below the plunger tip. Some of the sealing rings tested are, from left to right, top to bottom, a slotted Teflon ring with a spring insert, a pure Teflon o-ring, and a Teflon coated rubber o-ring.




Exhibit 8b. Disassembled test plunger. Screw-on tip was necessary because Teflon rings would not stretch over the plunger for mounting in a conventional groove.

Exhibit 8. Test plunger; identical in form to final plunger design.

DASH NO.	A	B	O-RING DASH NO. SEE NOTE 2
006	.235	.108	005
007	.266	.123	006
008	.297	.154	007
009	.329	.165	008
010	.360	.217	009
011	.422	.248	010
012	.485	.310	011
013	.550	.373	012
014	.613	.438	013
015	.675	.501	014
016	.738	.563	015
017	.800	.626	016
018	.863	.688	017
019	.925	.751	018
020	.993	.813	019
021	1.055	.881	020
022	1.118	.943	021
023	1.180	1.006	022
024	1.243	1.068	023
025	1.305	1.131	024
026	1.368	1.193	025
027	1.430	1.256	026
028	1.493	1.318	027
029	1.556	1.381	028
030	1.618	1.443	029
031	1.681	1.506	030
032	1.743	1.568	031
033	1.806	1.631	032
034	1.868	1.693	033
035	1.931	1.756	034
036	1.993	1.818	035
037	2.056	1.881	036
038	2.118	1.943	037
039	2.181	2.006	038
040	2.243	2.068	039
041	2.306	2.131	040
042	2.368	2.193	041
043	2.431	2.256	042
044	2.493	2.318	043
045	2.556	2.381	044
046	2.618	2.443	045
047	2.681	2.506	046
048	2.743	2.568	047
049	2.806	2.631	048
050	2.868	2.693	049
051	2.931	2.756	050
052	2.993	2.818	051
053	3.056	2.881	052
054	3.118	2.943	053
055	3.181	3.006	054
056	3.243	3.068	055
057	3.306	3.131	056
058	3.368	3.193	057
059	3.431	3.256	058
060	3.493	3.318	059
061	3.556	3.381	060
062	3.618	3.443	061
063	3.681	3.506	062
064	3.743	3.568	063
065	3.806	3.631	064
066	3.868	3.693	065
067	3.931	3.756	066
068	3.993	3.818	067
069	4.056	3.881	068
070	4.118	3.943	069
071	4.181	4.006	070
072	4.243	4.068	071
073	4.306	4.131	072
074	4.368	4.193	073
075	4.431	4.256	074
076	4.493	4.318	075
077	4.556	4.381	076
078	4.618	4.443	077
079	4.681	4.506	078
080	4.743	4.568	079
081	4.806	4.631	080
082	4.868	4.693	081
083	4.931	4.756	082
084	4.993	4.818	083
085	5.056	4.881	084
086	5.118	4.943	085
087	5.181	5.006	086
088	5.243	5.068	087
089	5.306	5.131	088
090	5.368	5.193	089
091	5.431	5.256	090
092	5.493	5.318	091
093	5.556	5.381	092
094	5.618	5.443	093
095	5.681	5.506	094
096	5.743	5.568	095
097	5.806	5.631	096
098	5.868	5.693	097
099	5.931	5.756	098
100	5.993	5.818	099
101	6.056	5.881	100
102	6.118	5.943	101
103	6.181	6.006	102
104	6.243	6.068	103
105	6.306	6.131	104
106	6.368	6.193	105
107	6.431	6.256	106
108	6.493	6.318	107
109	6.556	6.381	108
110	6.618	6.443	109
111	6.681	6.506	110
112	6.743	6.568	111
113	6.806	6.631	112
114	6.868	6.693	113
115	6.931	6.756	114
116	6.993	6.818	115
117	7.056	6.881	116
118	7.118	6.943	117
119	7.181	7.006	118
120	7.243	7.068	119
121	7.306	7.131	120
122	7.368	7.193	121
123	7.431	7.256	122
124	7.493	7.318	123
125	7.556	7.381	124
126	7.618	7.443	125
127	7.681	7.506	126
128	7.743	7.568	127
129	7.806	7.631	128
130	7.868	7.693	129
131	7.931	7.756	130
132	7.993	7.818	131
133	8.056	7.881	132
134	8.118	7.943	133
135	8.181	8.006	134
136	8.243	8.068	135
137	8.306	8.131	136
138	8.368	8.193	137
139	8.431	8.256	138
140	8.493	8.318	139
141	8.556	8.381	140
142	8.618	8.443	141
143	8.681	8.506	142
144	8.743	8.568	143
145	8.806	8.631	144
146	8.868	8.693	145
147	8.931	8.756	146
148	8.993	8.818	147
149	9.056	8.881	148
150	9.118	8.943	149
151	9.181	9.006	150
152	9.243	9.068	151
153	9.306	9.131	152
154	9.368	9.193	153
155	9.431	9.256	154
156	9.493	9.318	155
157	9.556	9.381	156
158	9.618	9.443	157
159	9.681	9.506	158
160	9.743	9.568	159
161	9.806	9.631	160
162	9.868	9.693	161
163	9.931	9.756	162
164	9.993	9.818	163
165	10.056	9.881	164
166	10.118	9.943	165
167	10.181	10.006	166
168	10.243	10.068	167
169	10.306	10.131	168
170	10.368	10.193	169
171	10.431	10.256	170
172	10.493	10.318	171
173	10.556	10.381	172
174	10.618	10.443	173
175	10.681	10.506	174
176	10.743	10.568	175
177	10.806	10.631	176
178	10.868	10.693	177
179	10.931	10.756	178
180	10.993	10.818	179
181	11.056	10.881	180
182	11.118	10.943	181
183	11.181	11.006	182
184	11.243	11.068	183
185	11.306	11.131	184
186	11.368	11.193	185
187	11.431	11.256	186
188	11.493	11.318	187
189	11.556	11.381	188
190	11.618	11.443	189
191	11.681	11.506	190
192	11.743	11.568	191
193	11.806	11.631	192
194	11.868	11.693	193
195	11.931	11.756	194
196	11.993	11.818	195
197	12.056	11.881	196
198	12.118	11.943	197
199	12.181	12.006	198
200	12.243	12.068	199

DASH NO.	G	R	E
006 TO 012	.075	.005/.015	.004
013 TO 028	.075	.005/.015	.005
111 TO 149	.112	.005/.015	.005
211 TO 243	.149	.010/.025	.006
316 TO 349	.221	.020/.030	.007
418 TO 437	.297	.020/.030	.010

NOTES:-

- DASH NUMBERS OF THIS Dwg. CORRESPOND TO DASH NUMBERS OF ARP 568 UNIFORM DASH NUMBERING SYSTEM FOR O-RINGS.
- O-RING DASH NUMBERS ARE SHOWN FOR SIZE DESIGNATION ONLY. SELECTION OF O-RING ELASTOMERIC COMPOUND SHOULD BE MADE TO SUIT THE FLUID AND TEMPERATURE ENVIRONMENT OF THE HYDRAULIC OR PNEUMATIC SERVICE IN WHICH THE SEAL WILL BE USED.
- THE SERIES "B" RING IS FURNISHED IN WEAR RESISTANT SPEC. 22-53 THE SERIES "B" RINGS ARE ALSO FURNISHED IN OTHER TFE MATERIALS UPON REQUEST.
- ALL DIMENSIONS SHOWN ARE IN INCHES.
- BORE DIAMETERS ARE PER MIL-P-5514 D.

Exhibit 9. Sealpruf Kap Seal Description.

MATERIAL, SIZE AND DESCRIPTION	
NO.	REV.
1	1

5-12068 D

NOTES:-

1. DASH NUMBERS OF THIS DWG. CORRESPOND TO DASH NUMBERS OF ARP 568 UNIFORM DASH NUMBERING SYSTEM FOR O-RINGS.
2. O-RING DASH NUMBERS ARE SHOWN FOR SIZE DESIGNATION ONLY. SELECTION OF O-RING ELASTOMERIC COMPOUND SHOULD BE MADE TO SUIT THE FLUID AND TEMPERATURE ENVIRONMENT OF THE HYDRAULIC OR PNEUMATIC SERVICE IN WHICH THE SEAL WILL BE USED.
3. THE SERIES "B" RING IS FURNISHED IN WEAR RESISTANT TURCON RINGS ARE (TYPE COMPOUND PER SPEC. ZZ-53) THE SERIES "B" ALSO FURNISHED IN OTHER TFE MATERIALS UPON REQUEST.
4. ALL DIMENSIONS SHOWN ARE IN INCHES.
5. ROBE TRANSFERS ARE PER MIL-B-5514 D.


Exhibit 9. Sealpruf Kap Seal Description.

Cap Seal Assembly

Spiral Wiper Ring

Tapered Edge for
Seal Entrance

Stud for Attaching
Cylinder Assembly
to Plunger Fork



Hexagon fitting

Metal Sleeve of
Cylinder

Joint of Metal
Sleeve to Glass
Cylinder

Plunger Tip

Cylinder
(Borosilicate
Precision Bore)

Plunger
(Delrin)

Exhibit 10a. Final
plunger-cylinder design.



Exhibit 10b. Final plunger and cylinder assembly.




Exhibit 10c. Disassembled plunger-cylinder unit. Seal used here is the Sealpruf Kap Seal.